

Portable Balances

CL Series

For Basic Weighing with a Simple 2-Button Keypad

- The CL is designed with superior RF protection to deliver stable results when used in production areas or near cell phones. Approvals for electrical and emissions include CE listed and FCC Class B Part 15 operation.
- The CL is designed with a large weighing platform and can adequately handle large samples including oversized letters.
- The CL has a low profile and storage is made easy with a stackable design. The integrated load cell lock protects the weighing cell from damage during storage or transportation.



Application Weighing

Display Liquid Crystal Display (LCD)

Operation AC Adapter (Not included) or 3 AA Batteries (Included)

Construction ABS Housing and Pan

Design Features Stacking design, low battery indicator, auto shut-off feature, mechanical and software overload protection, external push button calibration

Model	Capacity	Readability	Pan Size	Item No.	Price USD
CL201	200 g	0.1 g	Ø 4.75 in / Ø 120 mm	80010610	86
CL2000	2000 g	1 g	Ø 4.75 in / Ø 120 mm	80010611	86
CL5000	5000 g	1 g	Ø 4.75 in / Ø 120 mm	80010612	117
CL5000F	5000 g	1 g	Ø 4.75 in / Ø 120 mm	80010621	168

"F" model includes carrying case

Accessories	Item No.	Price USD
PWR,W,IN,6W,12VDC,S21x55x9x100,S	46001724	31
Pan Cover, SST, CL	80010623	11
PWR,ACCES,AC-PRONG	46001776	11