

Mechanical Scales & Balances

Dial-O-Gram[®] and Cent-O-Gram[®] 300 Series

High Precision and Remarkable Value

- Quick zeroing with zero adjust knob
- Vernier dial convenience on Dial-O-Gram (310 - 00)
- Self-aligning beam design

Display Easy to Read Deep Notched and Tiered Beams and Dial Plates

Construction All metal base and beam design with stainless steel pan

Design Features The three-point base, special floating bearing principle and beam design eliminates the need for level adjustment. Zeroing is quickly achieved by means of a zero adjust knob at the end of the beam. With aluminum pressure castings for the base and beam assembly, agate bearings, steel knife-edges, and stainless steel pan, everything about the Dial-O-Gram balance ensures a high standard of quality.



Model	Capacity	Readability	Beams	Pan Size	Pan Type	Item No.	Price USD
310-00	310 g	0.01 g	2 Beams with Dial	Ø 3.5 in / Ø 89 mm	Removable Pan	80000009	325
311-00	311 g	0.01 g	4 Beams	Ø 3.5 in / Ø 89 mm	Removable Pan	80000010	268

Accessories	Item No.	Price USD
Dust Cover, Dial-O-Gram Cent-O-Gram	80780005	26
Anti-Theft Device, SP PA AV NV	80850000	60

Harvard Trip[®]

Two Pan Design - Ideal for Comparative Weighing

- Built-in sliding masses
- Magnetic damping
- Floating agate bearings, counterbalancing knob for quick zeroing

Display Easy to Read Deep Notched and Tiered Beams

Construction Stainless steel plate

Design Features 2000g capacity with optional weight set. Built-in sliding masses for easy weighing, magnetic dampening to minimize oscillation and speed up weighing, self-aligning beam design and floating agate bearings, counterbalancing knob for quick zeroing, hardened steel edges



80000003



80000004

Model	Capacity	Readability	Weigh Beam Capacity	Beams	Pan Size	Item No.	Price USD
1450-SD	2000 g	0.1 g	10 g	Single Beam	Ø 6.2 in / Ø 152 mm	80000003	304
1550-SD	2000 g	0.1 g	210 g	2 Beams	Ø 6.2 in / Ø 152 mm	80000004	326
1560-SD	2000 g	0.1 g	210 g	2 Beams and Tare Beam	Ø 6.2 in / Ø 152 mm	80000005	340

Accessories	Item No.	Price USD	Accessories	Item No.	Price USD
Scoop, PP, 0.8l	80780006	26	Scoop and Counter Weight, 1.75l	80780016	39
Kit, Rod and Clamp	80780011	50	Anti-Theft Device, SP PA AV NV	80850000	60
Scoop and Counter Weight, Poly	80780014	37	Weight Set, 1kg-1g, ASTM 6, Cyl	80850115	369
Scoop and Counter Weight, SST	80780015	88	Dust Cover, Storage, Triple Beam	80850262	26