

Portable Balances

HH Series

Convenient, Compact And Portable Weighing Capability In The Palm Of Your Hand!

- Capacities from 120g to 320g are perfect for general weighing applications
- Backlit display for easy viewing in low lighting conditions and an energy-saving auto shut-off feature
- Durable ABS construction with protective snap on cover for storage and transport



Application	Weighing
Display	Liquid Crystal Display (LCD) with Backlight
Operation	2 AAA Batteries (Included)
Construction	ABS Housing, ABS Protective Cover, ABS Pan

Design Features Low battery indicator, overload and underload indicators, protective snap-on cover, external push button calibration, button activated luminescent display on HH120 and HH320

Model	Capacity	Readability	Pan Size	Item No.	Price USD
HH120	120 g	0.1 g	3.25 in x 3 in / 83 mm x 76 mm	71142842	166
HH120D	120 g	0.2 g	3.25 in x 3 in / 83 mm x 76 mm	71142841	139
HH320	320 g	0.1 g	3.25 in x 3 in / 83 mm x 76 mm	71142844	207

PS Series

Great Performance In A Unique Stylish Yet Functional Design!

- Sliding pan & cover protects display and keys and prevents load cell damage in transport
- Smaller size than most pocket scales allows a comfortable fit in your pocket, purse or briefcase.
- Fast and easy battery replacement with “pop-out” design



Application	Weighing
Display	Liquid Crystal Display (LCD)
Operation	2 Lithium Batteries (Included)
Construction	ABS Housing and Pan

Design Features Stability indicator, overload and underload indicators, low battery indicator, auto shut-off

Model	Capacity	Readability	Pan Size	Item No.	Price USD
PS121	120 g	0.1 g	2.5 in x 2.5 in / 64 mm x 64 mm	80104060	51
PS251	250 g	0.1 g	2.5 in x 2.5 in / 64 mm x 64 mm	80104061	67