

Moisture Analyzers

MB25

Redefining Value in Moisture Analysis!

- The MB25 features a bright backlit LCD and easy one-button setup for simple use and menu navigation.
- The MB25 utilizes a halogen heating element for swift drying of samples.
- The MB25's design makes it easy-to-clean and its compact size allows it to take up less space and fit into almost any work setting.



Application Weighing, Percent Moisture Determination, Percent Solid Determination

Display Liquid Crystal Display (LCD) with Backlight

Operation Power Cord (Included)

Communication RS232 Interface

Construction Halogen heat source, ABS housing, metal pan support, metal pan handler, replaceable in-use cover

Design Features 3 shut off criteria (manual, timed or auto), 1 heating profile (standard), 50° to 160° C heating range (5° C increments), displayed results (% moisture or % solids or weight, temperature, time)

Model	Capacity	Readability	Pan Size	Item No.	Price USD
MB25	110 g	5 mg	Ø 3.5 in / Ø 90 mm	83999803	3284

Accessories	Item No.	Price USD
Temp. Calibration Kit, MB Series	11113857	408
Holder, Pan, MB35 MB45	11113873	97
Printer, Impact, SF40A, AM	30064203	682
Printer, Thermal, STP-103, US	80251992	531
In-Use-Cover, MB23 MB25	80252475	25
Holder, Pan, MB Series	80252476	53
Cage, Sample, MB Series	80252477	69
Pan, Reusable, Set(3), 7mm, MB Series	80252478	57
Pan, Reusable, Set(3), 14mm, MB Series	80252479	60
Cable, ST103-AV DV EX MB PA TxxP	80252581	132
Cable, RS232, IBM 9P, AV DV MB Ranger	80500525	132
Cable, SF42-AV DV MB Ranger	80500571	132
Anti-Theft Device, MB AX EX	80850043	158
Software, Winwedge	80850080	374
Pan, Aluminum, Set(80), MB Series	80850086	42
Filter, Glass Fiber, Set(200)	80850087	70