

# Moisture Analyzers

## MB35

### High Performance Moisture Analysis!

- The MB35 has a dot matrix display that supports quick setup, shows percent moisture, actual weight, actual temperature and test time and also optimizes the viewing of results.
- MB35 includes a halogen heating lamp with an operating temperature range from 50° to 160° C in 5° increments for the quick drying of samples.
- The MB35 offers a straightforward one-step testing procedure that produces accurate results rapidly.



**Application** Weighing, Percent Moisture Determination

**Display** Dot Matrix Display with Backlight

**Operation** Power Cord (Included)

**Communication** RS232 with GLP/GMP Data Output

**Construction** Halogen heat source, Metal and ABS housing, metal pan support, metal pan handler, replaceable in-use cover

**Design Features** 3 shut off criteria (manual, timed or auto), 1 heating profile (standard), 50° to 160° C heating range (5° C increments), displayed results (% moisture or weight, temperature, time)

Model	Capacity	Readability	Pan Size	Item No.	Price USD
MB35	35 g	2 mg	Ø 3.5 in / Ø 90 mm	30079286	3824

Accessories	Item No.	Price USD
Temp. Calibration Kit, MB Series	11113857	408
Holder, Pan, MB35 MB45	11113873	97
Printer, Impact, SF40A, AM	30064203	682
Printer, Thermal, STP-103, US	80251992	531
Holder, Pan, MB Series	80252476	53
Cage, Sample, MB Series	80252477	69
Pan, Reusable, Set(3), 7mm, MB Series	80252478	57
Cable, ST103-AV DV EX MB PA TxxP	80252581	132
Cable, RS232, IBM 9P, AV DV MB Ranger	80500525	132
Cable, SF42-AV DV MB Ranger	80500571	132
Anti-Theft Device, MB AX EX	80850043	158
Software, Winwedge	80850080	374
In-Use-Cover, MB35 MB45	80850085	54
Pan, Aluminum, Set(80), MB Series	80850086	42
Filter, Glass Fiber, Set(200)	80850087	70