

Portable Balances

Navigator® NV

The Only Scale in its Class with Touchless Sensors That Free Up Your Hands!

- The Navigator Series is designed with two touchless sensors that will free up your hands so you can focus on handling samples.
- The Navigator Series is designed to accurately stabilize in less than one second to provide fast and accurate weighing results, improving operator efficiency, productivity and throughput.
- The Navigator Series is designed with a superior mechanical overload protection system that allows the scale to withstand loads up to 400% of its rated capacity.



- Application** Weighing, Percent Weighing, Parts Counting, Checkweighing, Display Hold, Accumulate/Totalize
- Display** Liquid Crystal Display (LCD) with Backlight
- Operation** AC Adapter (Included) or Internal Rechargeable Battery (Optional) or 4 C Batteries (Not included)
- Communication** Optional RS232 or USB or Ethernet Interface Kits with Integrated Cables
- Construction** ABS Housing, Stainless Steel Pan
- Design Features** Up front illuminated level indicator, removable stainless steel weighing platform, menu lockswitch, integral security bracket, integral shipping lock, stability indicator, overload and underload indicators, low battery indicator, auto shut-off

Model	Capacity	Readability	Pan Size	Item No.	Price USD
NV212	210 g	0.01 g	Ø 4.7 in / Ø 120 mm	83033081	313
NV211	210 g	0.1 g	Ø 4.7 in / Ø 120 mm	30135291	185
NV511	510 g	0.1 g	5.4 in x 7.5 in / 138 mm x 190 mm	83033083	214
NV1101	1100 g	0.1 g	5.4 in x 7.5 in / 138 mm x 190 mm	83033085	278
NV2101	2100 g	0.1 g	5.4 in x 7.5 in / 138 mm x 190 mm	83033087	362
NV4101	4100 g	0.2 g	5.4 in x 7.5 in / 138 mm x 190 mm	83033089	398
NV5101	5100 g	0.5 g	5.4 in x 7.5 in / 138 mm x 190 mm	83033091	438
NV4000	4000 g	1 g	5.4 in x 7.5 in / 138 mm x 190 mm	30135292	195

Accessories	Item No.	Price USD
Cable, PC, 9 Pin, NV AV	12106964	136
Auxiliary Display, TAD7, SP SPJ TA TAJ	80251395	45
Printer, Thermal, STP-103, US	80251992	531
Anti-Theft Device, SP PA AV NV	80850000	60
Software, Winwedge	80850080	374
RS232 Kit, SP TA Navigator	83032107	113
USB Kit, SPU TA NV NVL NVT	83032108	113
Ethernet Kit, NV NVL NVT	83032109	113
In-Use-Cover, NV	83032223	21
Carrying Case, NV	83032226	79