

# Portable Balances

## Navigator® XL

### The Only Scale in its Class with Touchless Sensors That Free Up Your Hands!

- The Navigator Series is designed with two touchless sensors that will free up your hands so you can focus on handling samples.
- The Navigator Series is designed to accurately stabilize in less than one second to provide fast and accurate weighing results, improving operator efficiency, productivity and throughput.
- The Navigator Series is designed with a superior mechanical overload protection system that allows the scale to withstand loads up to 400% of its rated capacity.



- Application** Weighing, Percent Weighing, Parts Counting, Checkweighing, Display Hold, Accumulate/Totalize
- Display** Liquid Crystal Display (LCD ) with Backlight
- Operation** AC Adapter (Included) or Internal Rechargeable Battery (Optional) or 4 C Batteries (Not included)
- Communication** Optional RS232 or USB or Ethernet Interface Kits with Integrated Cables
- Construction** ABS Housing, Stainless Steel Pan
- Design Features** Up front illuminated level indicator, removable stainless steel weighing platform, menu lockswitch, integral security bracket, integral shipping lock, stability indicator, overload and underload indicators, low battery indicator, auto shut-off

Model	Battery Power	Capacity	Readability	Pan Size	Item No.	Price USD
NVL511/1	•	510 g	0.1 g	7.6 in x 8 in / 194 mm x 203 mm	83033058	505
NVL1101/1	•	1100 g	0.1 g	7.6 in x 8 in / 194 mm x 203 mm	83033062	551
NVL2101/1	•	2100 g	0.1 g	7.6 in x 8 in / 194 mm x 203 mm	83033066	597
NVL5101/1	•	5100 g	0.5 g	7.6 in x 8 in / 194 mm x 203 mm	83033070	575
NVL10000/1	•	10000 g	1 g	7.6 in x 8 in / 194 mm x 203 mm	83033074	643
NVL20000/1	•	20000 g	1 g	7.6 in x 8 in / 194 mm x 203 mm	83033078	643

Accessories	Item No.	Price USD
Cable, PC, 9 Pin, NV AV	12106964	136
Auxiliary Display, TAD7, SP SPJ TA TAJ	80251395	45
Printer, Thermal, STP-103, US	80251992	531
Anti-Theft Device, SP PA AV NV	80850000	60
Software, Winwedge	80850080	374
Rechargeable Battery, NVL NVT	83032106	37
RS232 Kit, SP TA Navigator	83032107	113
USB Kit, SPU TA NV NVL NVT	83032108	113
Ethernet Kit, NV NVL NVT	83032109	113
In-Use-Cover, NVL	83032222	21
Carrying Case, NVL	83032225	79