

Analytical & Precision Balances

Pioneer® Analytical

The Economical Choice for Everyday Weighing!

- Pioneer's draftshield is constructed of all glass panels and three sliding doors, which are easy to remove for cleaning or storage.
- The up-front level indicator is located right next to the display and allows the user to quickly make sure the balance is level prior to each use.
- Pioneer is designed with selectable environmental settings and three filter modes in order to ensure accurate operation in a variety of settings.



Application Weighing, Parts Counting, Percentage Weighing, Density Determination and Dynamic Weighing Applications

Display Backlight Liquid Crystal Display (LCD)

Operation AC Adapter (Included)

Communication RS232 (included)

Construction Solid diecast metal base, stainless steel pan

| Model | InCal™ | Capacity | Readability | Pan Size | Item No. | Price USD |
|--------|--------|----------|-------------|--------------------|----------|-----------|
| PA84 | | 85 g | 0.1 mg | Ø 3.5 in / Ø 90 mm | 30208440 | 1605 |
| PA84C | • | 85 g | 0.1 mg | Ø 3.5 in / Ø 90 mm | 30208454 | 1805 |
| PA124 | | 120 g | 0.1 mg | Ø 3.5 in / Ø 90 mm | 30208441 | 1926 |
| PA124C | • | 120 g | 0.1 mg | Ø 3.5 in / Ø 90 mm | 30208455 | 2126 |
| PA224 | | 220 g | 0.1 mg | Ø 3.5 in / Ø 90 mm | 30208442 | 2140 |
| PA224C | • | 220 g | 0.1 mg | Ø 3.5 in / Ø 90 mm | 30208456 | 2340 |

| Accessories | Item No. | Price USD |
|--|----------|-----------|
| Cable, RS9 M/F, 1.5m | 12122603 | 132 |
| Printer, Impact, SF40A, AM | 30064203 | 682 |
| Auxiliary Display, PAD7, PA AX EX | 80251396 | 45 |
| Printer, Thermal, STP-103, US | 80251992 | 531 |
| Cable, RS232, CBM910-AV DV EX MB PA TxxP | 80252571 | 132 |
| Cable, ST103-AV DV EX MB PA TxxP | 80252581 | 132 |
| Density Kit, Solids | 80253384 | 533 |
| Cable, 25 Pin-9 Pin, PC-TxxP | 80500524 | 132 |
| Cable, RS232, IBM 9P, AV DV MB Ranger | 80500525 | 132 |
| Adapter, 25 Pin-9 Pin, SF42 SRP275 | 80500572 | 45 |
| Anti-Theft Device, SP PA AV NV | 80850000 | 60 |
| Software, Winwedge | 80850080 | 374 |
| In-Use-Cover, PA | 83020221 | 31 |