

Compact Bench Scales

Ranger® 3000

The Best Compact Scale In Its Class!

- Ranger 3000 produces accurate results within one second. The rapid stabilization virtually eliminates your wait time and allows you to move at a pace determined by you, not your weighing equipment.
- Ranger 3000 has the largest display in its class and is equipped with Checkweighing LED's and backlit LCD displays which makes it functional in almost any working condition.
- The most user friendly scale on the market with easy setup ensures that your valuable time is not wasted and you're working as proficiently as possible.



Application Weighing, Parts Counting, Checkweighing, Percent Weighing, Dynamic Weighing, Display Hold, Accumulation

Display Liquid Crystal Display (LCD) with Backlight

Operation AC Power (Included) or Rechargeable Battery (Included)

Communication Easy access communication ports including a standard RS232 and an optional second port of either RS232, USB or Ethernet

Construction ABS Housing, Stainless Steel Platform

Design Features 210 hours of battery life, 3 LED (yellow, green, red) Checkweighing lights with selectable operation and audible signal settings, integral weigh-below hook, transportation handle, sealed front panel, menu lockswitch, up front level indicator, adjustable leveling feet, selectable environmental and auto-print settings, stability indicator, overload and underload indicators, low battery indicator, auto shut-off, auto tare

Model	Capacity	Readability	Certified Readability	Platform Size	Item No.	Price USD
R31P1502	3 lb / 1.5 kg	0.0001 lb / 0.05 g	0.001 lb / 0.5 g	8.9 in x 11.8 in / 225 mm x 300 mm	30031707	504
R31P3	6 lb / 3 kg	0.0002 lb / 0.1 g	0.002 lb / 1 g	8.9 in x 11.8 in / 225 mm x 300 mm	30031708	504
R31P6	15 lb / 6 kg	0.0005 lb / 0.2 g	0.005 lb / 2 g	8.9 in x 11.8 in / 225 mm x 300 mm	30031709	504
R31P15	30 lb / 15 kg	0.001 lb / 0.5 g	0.01 lb / 5 g	8.9 in x 11.8 in / 225 mm x 300 mm	30031710	504
R31P30	60 lb / 30 kg	0.002 lb / 1 g	0.02 lb / 10 g	8.9 in x 11.8 in / 225 mm x 300 mm	30031711	504

Accessories	Item No.	Price USD
Ethernet Kit, R31 RC31 R71 V71	30037447	104
RS232 Kit, R31 RC31 R71 V71	30037448	104
USB Kit, R31 RC31 V71	30037449	104
In-Use-Cover, R21 RC21 R31 RC31	30037451	24
Printer, Impact, SF40A, AM	30064203	682
Printer, Thermal, STP-103, US	80251992	531
Cable, RS232, IBM 9P, AV DV MB Ranger	80500525	132
Carrying Case, TR TC R31 TC31 V71	80850084	259