

Printers

SF40A

OHAUS Value in a Full-Featured, Portable Printer

- The statistics and totalization modes and real time clock complement OHAUS weighing instruments' application modes and offer functionality that assists the weight determination and reporting processes.
- SF40A helps support the efficiency of the weighing and recording processes. RS232 interface and PC configuration software make setting up and using the printer fast and simple.
- With the ability to print in 13 different languages, SF40A can print in the language of the connected device to support universal operation.



- Application** Statistics, Totalization
- Display** Not Applicable
- Operation** AC Adapter (Included)
- Communication** Easy Access RS232 Communication Port
- Construction** ABS Housing
- Design Features** Status indicator, 40 characters per line, compatible with continuous labelpaper, perform zero and tare operations directly from the printer

Model	Item No.	Price USD
SF40A	30064203	682

STP103

Economical Palm-Sized Thermal Printer

- Ultra compact and features low power consumption
- 42 column printing at a speed of approximately 11.8 lines per second
- Simple paper loading and suspends printing when printing paper becomes low



- Application** Not Applicable
- Display** Not Applicable
- Operation** AC Adapter (Included)
- Communication** Easy Access RS232 Communication Port
- Construction** ABS Housing
- Design Features** Status indicator, paper feed, graphics printing capability

Model	Item No.	Price USD
Printer, Thermal,STP103, 120vac US	80251992	531