

Mechanical Scales & Balances

School Balance

Unleash Your Young Scientists!

- Large, washable weighing pans
- 18 piece stackable plastic mass set
- Dual pointers to allow viewing results from both sides



Display Plastic pointer indicates balanced condition

Construction ABS Base, Beam and Weighing Pans

Model	Capacity	Readability	Pan Size	Item No.	Price USD
SB1200	2000 g	0.5 g	5.5 in x 5.5 in / 140 mm x 140 mm	80000029	50
SB1200 Measurement Bundle	2000 g	0.5 g	5.5 in x 5.5 in / 140 mm x 140 mm	80252137	337

Primer[®] Balance

Unleash Your Young Scientists!

- Activity guide with 16 exercises
- Easy to use zero adjustment
- Dual pointers to allow viewing results from both sides



Display Plastic pointer indicates balanced condition

Construction ABS Base, Beam and Weighing Buckets

Model	Capacity	Readability	Pan Size	Item No.	Price USD
80410-00	2000 g	1 g	4 in x 4 in / 102 mm x 102 mm	80000024	36