

Portable Balances

Scout® Pro

Hands Down The Industry's Best Portable Balance!

- Mechanical and software protection provides superior resistance to shock and overloading and prevents damage to the weighing cell.
- Easily connect to peripheral devices through the installable RS232 or USB interface.
- Scout Pro offers Enhanced Parts Counting, Percent Weighing, Totalization, and Display Hold.



Application Weighing, Parts Counting with Average Piece Weight Recalculation, Percentage Weighing, Totalization, Display Hold

Display Liquid Crystal Display (LCD)

Operation AC Adapter (Included) or 4 AA Batteries (Not included)

Communication Optional RS232 or USB Interface with integrated cable

Construction ABS Housing, Stainless Steel Pan

Design Features Software lockswitch, integral security bracket, integral shipping lock, integral weigh below hook, sealed front panel and molded spill ring, stability indicator, overload and underload indicators, low battery indicator, auto shut-off, span calibration mass included on 200g, 400g and 600g models

Model	Capacity	Readability	Pan Size	Item No.	Price USD
SP202	200 g	0.01 g	Ø 4.7 in / Ø 120 mm	30064412	428
SP402	400 g	0.01 g	Ø 4.7 in / Ø 120 mm	30064413	600
SP602	600 g	0.01 g	Ø 4.7 in / Ø 120 mm	30064419	702
SP401	400 g	0.1 g	Ø 4.7 in / Ø 120 mm	30064414	236
SP601	600 g	0.1 g	5.6 in x 6.5 in / 142 mm x 165 mm	30064415	349
SP2001	2000 g	0.1 g	5.6 in x 6.5 in / 142 mm x 165 mm	30064416	441
SP4001	4000 g	0.1 g	5.6 in x 6.5 in / 142 mm x 165 mm	30064417	646
SP6001	6000 g	0.1 g	5.6 in x 6.5 in / 142 mm x 165 mm	30064420	728
SP6000	6000 g	1 g	5.6 in x 6.5 in / 142 mm x 165 mm	30064418	366

Accessories	Item No.	Price USD
Printer, Impact, SF40A, AM	30064203	682
In-Use-Cover, Set(3), SP	71168909	37
Printer, Thermal, STP-103, US	80251992	531
Cable, RS232,CBM910-AV DV EX MB PA TxxP	80252571	132
Cable, ST103-AV DV EX MB PA TxxP	80252581	132
Adapter, 25 Pin-9 Pin, SF42 SRP275	80500572	45
Anti-Theft Device, SP PA AV NV	80850000	60
Software, Winwedge	80850080	374
RS232 Kit, SP TA Navigator	83032107	113
USB Kit, SP TA NV NVL NVT	83032108	113
Density Kit, SP, SPJ	71147378	13