

Mechanical Scales & Balances

Triple Beam 700 Series

The Standard by Which All Others Are Measured!

- Superior magnetic damping provides stability reading three times faster than any other three-beam balance on the market
- Choice of 5 different weighing pans to fit any application
- Counterbalancing knob for quick zeroing, self-aligning beam design

Display Easy to Read Deep Notched and Tiered Beams

Construction All metal base and beam design with stainless steel pan

Design Features Positive poise positioning for repeatable results, self-aligning beam design, magnetic damping to minimise oscillation and speed weighing



80000012



80000034

Model	Capacity	Readability	Pan Size	Pan Height	Pan Type	Item No.	Price USD
720-00	610 g	0.1 g	12 in x 6 in / 30.5 mm x 152 mm	2.75 in/70 mm	Polypropylene Scoop	80000032	264
720-50	610 g	0.1 g	12 in x 7.4 in / 305 mm x 188 mm	1.7 in/43 mm	Stainless Steel Scoop	80000033	306
710-00	610 g	0.1 g	Ø 6 in / Ø 152 mm	—	Removable Pan	80000011	223
750-50	610 g	0.1 g	Ø 6 in / Ø 152 mm	—	Fixed Pan	80000012	163
760-00	610 g	0.1 g	Ø 6 in / Ø 152 mm	—	Fixed Pan	80000013	231
710-T0	610 g	0.1 g	Ø 6 in / Ø 152 mm	—	Removable Pan	80000031	273
750-50 with Virtual Lab	610 g	0.1 g	Ø 6 in / Ø 152 mm	—	Fixed Pan	80104165	210
730-00	610 g	0.1 g	Ø 9 in / Ø 229 mm	6 in/152 mm	Animal Container	80000034	330

Accessories	Item No.	Price USD
Subject Box	80000053	70
Scoop and Counter Weight	80780001	41
Scoop, PP, 0.8l	80780006	26
Kit, Rod and Clamp	80780011	50
Scoop and Counter Weight, Poly	80780014	37
Scoop and Counter Weight, SST	80780015	88
Scoop and Counter Weight, 1.75l	80780016	39
Dust Cover, Triple Beam	80780017	26
Weight Set, 2x1kg, 1x500g	80780108	52
Anti-Theft Device, SP PA AV NV	80850000	60
Dust Cover, Storage, Triple Beam	80850262	26
Clamp, Gas Cylinder	80850312	67



80000011



80000032



80000033



80104165