

# Mechanical Scales & Balances

## Triple Beam Junior

### OHAUS Quality at an Affordable Price

- Notched and tiered beams for easy reading
- Security Features (lockable connections)
- Zero adjust compensator ensures correct zero

**Display** Easy to Read Deep Notched and Tiered Beams

**Construction** All metal base and beam design with stainless steel pan

**Design Features** Notched and tiered beams for easy reading, magnetic dampening and zero adjust compensator ensures correct zero.



Model	Capacity	Readability	Pan Size	Item No.	Price USD
TJ611	610 g	0.1 g	Ø 5.9 in / Ø 150 mm	80252052	147

Accessories	Item No.	Price USD
Dust Cover, Triple Beam	80780017	26
Anti-Theft Device, SP PA AV NV	80850000	60
Dust Cover, Storage, Triple Beam	80850262	26

## Triple Beam Dial-O-Gram® 1600 Series

### The Standard by Which All Others Are Measured!

- Superior magnetic damping provides three times faster stability reading than any other three-beam balance on the market
- Vernier dial convenience
- Counterbalancing knob for quick zeroing, self-aligning beam design

**Display** Easy to Read Deep Notched and Tiered Beams and Dial Plates

**Construction** Metal Base and Beam, Stainless Steel Pan

**Design Features** Counterbalancing knob for quick zeroing, positive poise positioning for repeatable results, self-aligning beam design, magnetic dampening to minimise oscillation and speed weighing.



80000037



80000006

Model	Capacity	Readability	Pan Size	Pan Type	Item No.	Price USD
1650-00	610 g	0.1 g	Ø 6 in / Ø 152 mm	Fixed Pan	80000006	294
1610-00	610 g	0.1 g	Ø 6 in / Ø 152 mm	Removable Pan	80000036	338
1650-W0	2610 g	0.1 g	Ø 6 in / Ø 152 mm	Fixed Pan	80000037	338

Accessories	Item No.	Price USD
Scoop, Flat Bottomed, SST, 0.6l	80250400	50
Scoop and Counter Weight	80780001	41
Scoop, PP, 0.8l	80780006	26
Kit, Rod and Clamp	80780011	50
Scoop and Counter Weight, Poly	80780014	37
Scoop and Counter Weight, SST	80780015	88
Scoop and Counter Weight, 1.75l	80780016	39
Dust Cover, Triple Beam	80780017	26

Accessories	Item No.	Price USD
Weight Set, 2x1kg, 1x500g	80780108	52
Anti-Theft Device, SP PA AV NV	80850000	60
Dust Cover, Storage, Triple Beam	80850262	26
Clamp, Gas Cylinder	80850312	67