

Compact Bench Scales

Valor® 3000

Rugged, Full Stainless Steel Design In A High Performance, Easy-To-Use, Portable Scale.

- High resolution models available!
- NEMA 4x/IP65 models available!
- NSF Certified, USDA-AMS Accepted, supports HACCP-Certified systems



- Application** Weighing, Percent Weighing, Specific Gravity
- Display** Liquid Crystal Display (LCD) with white LED backlight
- Operation** AC Adapter (Included) or 4 C (LR14) Batteries (Not Included)
- Construction** Stainless steel housing and stainless steel pan
- Design Features** Adjustable leveling feet and bubble on XW and N models, load cell lock, user adjustable filtering and optional food accessories



Model	Protection Rating	Capacity	Readability	Certified Readability	Platform Size	Item No.	Price USD
V31XH202	Not Applicable	0.44 lb / 200 g	0.00005 lb / 0.01 g	—	Ø 4.7 in / Ø 120 mm	83998130	592
V31XH402	Not Applicable	0.88 lb / 400 g	0.00005 lb / 0.01 g	—	Ø 4.7 in / Ø 120 mm	83998178	766
V31XW301	IP65; NEMA 4X	0.66 lb / 300 g	0.0005 lb / 0.1 g	—	6.2 in x 5.7 in / 158 mm x 146 mm	83998134	479
V31X501	Not Applicable	1.1 lb / 500 g	0.0005 lb / 0.1 g	—	Ø 4.7 in / Ø 120 mm	83998176	446
V31XH2	Not Applicable	4.4 lb / 2 kg	0.0005 lb / 0.1 g	—	6.2 in x 5.7 in / 158 mm x 146 mm	83998131	592
V31XH4	Not Applicable	8.8 lb / 4 kg	0.0005 lb / 0.1 g	—	6.2 in x 5.7 in / 158 mm x 146 mm	83998177	813
V31X3	Not Applicable	6.6 lb / 3 kg	0.005 lb / 1 g	—	6.2 in x 5.7 in / 158 mm x 146 mm	83998132	405
V31X6	Not Applicable	13.2 lb / 6 kg	0.005 lb / 1 g	—	6.2 in x 5.7 in / 158 mm x 146 mm	83998133	405

Model	Protection Rating	Capacity	Readability	Certified Readability	Platform Size	Item No.	Price USD
V31X3N	Not Applicable	6.6 lb / 3 kg	0.005 lb / 0.5 g	0.005 lb / 1 g	6.2 in x 5.7 in / 158 mm x 146 mm	83999992	473
V31XW3	IP65; NEMA 4X	6.6 lb / 3 kg	0.005 lb / 0.5 g	0.005 lb / 1 g	6.2 in x 5.7 in / 158 mm x 146 mm	83998135	479
V31XW6	IP65; NEMA 4X	13.2 lb / 6 kg	0.005 lb / 1 g	0.005 lb / 2 g	6.2 in x 5.7 in / 158 mm x 146 mm	83998136	479
V31X6N	Not Applicable	13.2 lb / 6 kg	0.005 lb / 1 g	0.005 lb / 2 g	6.2 in x 5.7 in / 158 mm x 146 mm	83999996	473



Accessories	Item No.	Price USD
Pan, SST, Sandwich, V31	80251147	55
Holder, SST, French Fries, V31	80251148	55
Holder, SST, Ice Cream Cone, V31	80251149	55
Carrying Case, V31	80251216	88